

Get better with ActivaScrew™ Interference

High Usability for Various Indications

ActivaScrew™ Interference offers

- HIGH USABILITY FOR VARIOUS INDICATIONS
- WIDE SELECTION OF DIFFERENT SIZES
- EASY INSERTION ALONG THE GUIDE WIRE (K-WIRE)
- EASY POSITION IMAGING DURING THE OPERATION
- CENTER DRIVEN TORX DRIVER FOR OPTIMAL INSERTION
- HIGH STRENGTH PROPERTIES OFFER EASY INSERTION AND SAFE MEDICAL USE
- SELF-LOCKING SL™
- ISOELASTICITY
- NO MR IMAGING DISTURBANCE, MR SAFE
- ENABLES DRILLING THROUGH THE PRODUCT IN NEED OF RE-OPERATION
- NO NEED FOR REMOVAL OPERATION
- IMPLANT IS SUPPLIED GAMMA STERILIZED



The ActivaScrew™ Interference bioabsorbable screws are constructed of bioabsorbable poly(L-lactide-co-glycolide) (PLGA).

ActivaScrew™ Interference Instruments

Product Code Description

B-II-9031	Torx driver for 7/8/9/10 mm interference screw
B-II-8010	Notcher for 7/8/9/10 mm interference screw
B-II-1624	K-wire Ø1.6 mm, length 240 mm, for 7/8/9/10 mm interference screw
B-II-1624B	K-wire Ø1.6 mm, length 240 mm, blunt ends, for 7/8/9/10 mm interference screw



Manufacturer
Bioretec Ltd.
Hermiankatu 22, Modulight Building
FI-33720 Tampere, Finland
Tel. +358 20 778 9500
Fax +358 3 317 0225
sales@bioretec.com
www.bioretec.com

Distributor

Surgical Technique For ACL Reconstruction

1 Drilling of Tibial and Femoral Tunnels

After harvesting and/or preparation of the soft tissue graft, prepare the tibial and femoral tunnels by drilling. Use an appropriate drill bit according to the graft size. Use a sufficient drill guide.



2 Graft Insertion

Insert the graft into the drilled tunnels. If a soft tissue graft is used, both ends should be whip stitched 3 cm to prevent twisting or wrapping during screw insertion.



3 Notching

Insert the proper size guide wire into the femoral tunnel and use the cannulated Bioretac Notcher (B-II-8010) to prepare the insertion point of screw entrance.



4 Femoral Interference Fixation

Using Bioretac Interference Torx driver (B-II-9031) insert the ActivaScrew™ Interference along the guide wire (B-II-1624 or B-II-1624B) into the femoral tunnel.

Typically, 7 - 10 mm screws are used for ACL reconstruction.



5 Tibial Interference Fixation

Tension the graft from the tibial side and insert the tibial screw using Bioretac Interference Torx driver (B-II-9031). If needed, Bioretac Notcher (B-II-8010) can be used also with the tibial tunnel.

