


GET BETTER
with ActivaPin™ you will



Bioretec products have their roots in research on biomaterials and implant technology started by the team led by Professor Pertti Törmälä at the Tampere University of Technology in 1977, first in the world.

So far, millions of patients have benefited of earlier innovations of our key people. **Bioretec's** products are developed under customer guidance to fulfill the highest demands of safety, user-friendliness and cost-efficiency.

Bioretec is devoted to supply solutions, which satisfy surgeon requirements and benefit patients. Financially feasible and responsible materials are the top priority and an essential part of the high-tech solutions we offer.

Implants that do more – It's Bioretec

www.bioretec.com



Manufacturer

Bioretec Ltd.
Hermiankatu 22, Modulight Building
FI-33720 Tampere, Finland
Tel. +358 20 778 9500
Fax +358 3 317 0225
sales@bioretec.com
www.bioretec.com

Distributor



B-L-APB1, Rev.1.0

www.bioretec.com

Properties and techniques of the ActivaPin™

bioretec

GET BETTER with ActivaPin™

The new generation of bioabsorbable implants




Figure 1: Cross section of the ActivaPin™ 1.5 mm (left), indicating maximum and core diameter measurements. The illustration on the right demonstrates the relationship of the device to a 1.5 mm drill hole (right).

Product Range

ActivaPin™

REMARK

Prior to using Bioretec Implants examine thoroughly the Instructions For Use - inside each product package.

ORDERING INFORMATION

For more information please contact your nearest distributor or our customer service at sales@bioretec.com. Valid product range available from our website: www.bioretec.com.

ActivaPin™

Product Code	Description
B-AP-1520	1.5 mm x 20 mm
B-AP-1530	1.5 mm x 30 mm
B-AP-1540	1.5 mm x 40 mm
B-AP-1550	1.5 mm x 50 mm
B-AP-1560	1.5 mm x 60 mm
B-AP-1570	1.5 mm x 70 mm
B-AP-2020	2.0 mm x 20 mm
B-AP-2030	2.0 mm x 30 mm
B-AP-2040	2.0 mm x 40 mm
B-AP-2050	2.0 mm x 50 mm
B-AP-2060	2.0 mm x 60 mm
B-AP-2070	2.0 mm x 70 mm
B-AP-2750	2.7 mm x 50 mm
B-AP-2770	2.7 mm x 70 mm
B-AP-3270	3.2 mm x 70 mm

ActivaPin™ Instruments

Product Code	Description
B-IP-1500	Applicator for 1.5 mm ActivaPin™
B-IP-2000	Applicator for 2.0 mm ActivaPin™
B-IP-2700	Applicator for 2.7 mm ActivaPin™
B-IP-3200	Applicator for 3.2 mm ActivaPin™
B-IP-1503	Drill Bit Ø 1.5 mm (85 mm) with quick coupling
B-IP-2003	Drill Bit Ø 2.0 mm (100 mm) with quick coupling
B-IP-2703	Drill Bit Ø 2.7 mm (100 mm) with quick coupling
B-IS-4510	Drill Bit Ø 3.2 mm (120 mm) with quick coupling
B-IP-1501	K-wires (10 pcs) for 1.5 mm ActivaPin™
B-IP-2001	K-wires (10 pcs) for 2.0 mm ActivaPin™
B-IP-2701	K-wires (10 pcs) for 2.7 mm ActivaPin™
B-IP-3201	K-wires (10 pcs) for 3.2 mm ActivaPin™

Accessories

B-HTC-1000	High Temperature Cautery
B-IP-4003	Sterilization tray

ActivaPin™

Patents: EP 1864616, EP 1902680, FI 124190. EP, US, FI Patent Pending